# MICRO STONE DISC MILL

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销售  
Sales


售后  
After sales

## 设备原理 Equipment principle

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The upper plate is fixed and the lower plate rotates at a high speed by using the ceramic grinding plate of the upper and lower groups. The resulting kneading, shearing, depolymerization and dispersion achieve the crushing effect.

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Equipment name: Micro stone disc mill

Use: For superfine crushing of raw materials

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
【Application】

Chemicals: Battery material Lithium hydroxide Graphite Nickel oxide Metal oxide etc.

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Internal cavity: White alumina as a whole

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Scale instrument imported from Japan

Precision adjustable: 1 ~ 100 small grids

M1 calculation scale:  $0.25/100 = 0.0025$

$0.0025 \times 4 = 0.01\text{mm}$   $0.01\text{mm} = 10\mu\text{m}$

Where contact noise was heard between the fixed and rotary grinders.

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A jacket is set in the grinding chamber and the upper cover, and cold water is circulated to reduce the grinding temperature.

型号 Model	XRS-150	XRS-250	XRS-360
磨盘转速 Disk speed	1,800 r/min	3,000 r/min	1,500 r/min
电动机 Electric motor	3.7kW	11kW	22kW
磨盘尺寸 Disk size	$\phi 150\text{mm}$	$\phi 250\text{mm}$	$\phi 360\text{mm}$
机器尺寸 Equipment dimensions	W560×L350×H950mm	W890×L400×H1,300mm	W700×L960×H1,800mm
处理能力 Processing capacity	30 ~ 100kg	100 ~ 1,000kg	300 ~ 2,500kg
设备重量 Equipment weight	120kg	250kg	400kg

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Enterprises all over the world can manufacture stone plate mills  
Only our company has full ceramic lining!  
Start to monopolize the battery industry!

- ❖所有原材料剪断·磨碎·微粉化·分散·乳化
- ❖粉碎物料粒径是通过调节磨盘的间距来获取
- ❖通过强力撞击·离心力·剪力进行粉碎
- ❖即使在运行中也可调整磨盘间距
- ❖Cutting, grinding, micronization, dispersion and emulsification of all raw materials.
- ❖The particle size of crushed materials is obtained by adjusting the spacing of grinding discs.
- ❖Crushing by strong impact, centrifugal force and shearing force.
- ❖The grinding disc spacing can be adjusted even during operation.

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4 Chinese patent authorization  
Announcement No.

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Crushing is done by the equipment alone.  
No need for cyclones, bag filters, blowers, etc.  
Highly effective in reducing capital investment.

**其他公司不具备的优点和特长**  
Advantages that other companies do not have

- ❖易清洗、拆卸简单,维修简单
- ❖设计新颖,操作方便
- ❖噪音低、破碎率高
- ❖干粉原料·湿粉原料两用
- ❖粉碎腔体及上盖内设有水冷套
- ❖可用于实验室
- ❖即使少量原料做实验可100%回收
- ❖符合CE·KCs·UL·CCC安全标准的设计
- ❖Easy to clean, easy to disassemble, easy to maintain.
- ❖New design, easy to operate.
- ❖Low noise and high crushing rate.
- ❖Dry powder raw material and wet powder raw material are both used.
- ❖Water cooling jacket in the crushing chamber and upper cover.
- ❖Can be used in laboratory.
- ❖Even a small amount of raw material can be 100% recycled for experiments.
- ❖Designed to conform CE/KCs/UL/CCC safety standards.

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The control box adopts a frequency conversion system,  
the number of rotations is adjustable,  
and the temperature is considerable.

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
Easy to disassemble triangular wing, easy to clean.

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Less stagnant materials in the cavity.

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
The air flow drives the wings to promote  
the discharge of raw materials.

# 异物混入防止对策 Measures to prevent contamination


① 漏斗内原料接粉部  
特氟龙涂层  
① Raw material powder receiving part  
in hopper teflon coating

③ 内衬氧化铝  
③ Lined with alumina

④ 两个桨叶可选用耐磨性树脂  
或碳化钨制作  
④ Two blades are made of wear-resistant resin  
or tungsten carbide

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
⑤ Fixing screw (all tungsten carbide alloy)  
Fabrication or alumina ceramics  
Alumina purity ≥ 99.9%  
Very wear resistant

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
Problem Measures to prevent damage to pulverizers due to contamination of raw materials with metallic foreign matter.

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### ◆ Working principle

- ◆ Installed on the top of the micro grinder to protect the ceramic parts from foreign matters.
- ◆ Fine powder is supplied to the micro pulverizer through the vibration of the air compressor vibrator.
- ◆ Foreign matters such as metal parts remain on the screen.
- ◆ Screen mesh can be freely selected.
- ◆ Raw material supply port and outlet flange can also choose the required size.
- ◆ Removal of foreign matters is manual.

即便是1颗螺钉混入原料  
对微粉碎机的磨盘是致命的打击!  
如果从以前开始有这样的机器...  
而且比较振动筛压倒性的小型!  
If screws were mixed in with the raw material,  
the rotor of the crusher would break into pieces.  
If only we had such a machine before...  
And it's much smaller than a vibrating sieve!

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


Patent No. : CN108452903A  
Our company's high-density, non porous stone mortar  
Excellent sharpness and handling ability to prevent  
bacteria from entering and breeding.

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


April 2022 XRS-360 Type  
Achieved D50=5.9μm±0.3μm!  
Opens a new era of stone disk grinding!  
Production capacity of 700kg/hr!

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中国科学院苏州纳米研究所产学研基地  
Suzhou Institute of Nanotechnology, Chinese Academy of Sciences, Designation of Joint Research Centers


中南大学合作设立博士生工作站  
Collaboration with Central South University to establish a doctoral internship

ISO9001认证  
ISO9001 certification

ISO14001认证  
ISO14001 certification

专利证书  
Patent certificate

欧洲CE认证  
CE certification in Europe

韩国KCs认证  
KCs certification in Korea

进出口业务许可证  
License of import and export business

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防爆设备安装认证  
Explosion-proof equipment installation certification

防爆设备维护认证  
Explosion-proof equipment maintenance certification

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XiRan XiRan Technology (Jiangsu) Co., Ltd.